

IHC Double staining kit(Polymer-HRP)(DAB/AEC, For Mouse/Rabbit Primary Antibodies)

Catalog No.: PA10002

Basic Information

Product name	IHC Double staining kit(Polymer-HRP)(DAB/AEC, For Mouse/Rabbit Primary Antibodies)
Sizes	6 mL, 30 mL, 100 mL
Storage	2-8°C
Shipping	Shipped with ice pack
Validity	12 months

Product Introduction

EnkiLife's dual immunohistochemistry detection system enables simultaneous detection of two distinct antigens on the same tissue section. The Polymer-HRP-based IHC detection system employs advanced polymer technology to conjugate horseradish peroxidase (HRP) with specific antibodies, forming a stable polymer complex. This unique labeling approach significantly enhances signal intensity and detection sensitivity while minimizing background staining, offering user-friendly operation and robust stability. The kit includes Polymer-HRP-labeled anti-mouse and anti-rabbit secondary IgGs, paired with DAB (brown) and AEC (red) chromogenic substrates, to achieve dual-antigen localization and visualization.

Product Components

Components	6 mL	30 mL	100 mL
Reagent 1: 3 % H ₂ O ₂	6 mL	30 mL	100 mL
Reagent 2: Goat serum blocking solution	12 mL	60 mL	200 mL
Reagent 3: Polymer-HRP-conjugated goat anti-mouse IgG	6 mL	30 mL	100 mL
Reagent 4: Polymer-HRP-conjugated goat anti-rabbit IgG	6 mL	30 mL	100 mL
Reagent 5: DAB concentrate (50×)	120 µL	600 µL	2 mL
Reagent 6: DAB diluent	6 mL	30 mL	100 mL
Reagent 7: AEC chromogen A	300 µL	1.5 mL	5 mL

IHC Double staining kit(Polymer-HRP)(DAB/AEC, For Mouse/Rabbit Primary Antibodies)

Catalog No.: PA10002

Reagent 8: AEC chromogen B	300 µL	1.5 mL	5 mL
Reagent 9: AEC chromogen C	6 mL	30 mL	100 mL

Experimental procedure

1. Deparaffinization and rehydration: Dewax formalin-fixed paraffin-embedded (FFPE) sections in xylene and rehydrate through graded ethanol.
2. Antigen retrieval: Perform heat-induced epitope retrieval (HIER) or enzymatic retrieval as required by the primary antibodies.
3. Endogenous peroxidase quenching: Incubate slides with 3 % H₂O₂ (Reagent 1) for 15–20 min at RT.
4. Wash: 3 × 3–5 min with PBST.
5. Blocking: Cover tissue with goat serum blocking solution (Reagent 2) and incubate 30 min at 37 °C.
6. Primary antibody incubation: After determining the optimal dilution, evenly overlay the tissue section with the working antibody solution and incubate overnight at 4 °C or for 1 h at 37 °C.
7. Wash: 3 × 3–5 min PBST.
8. Secondary antibody incubation: Add Polymer-HRP-labeled goat anti-mouse IgG (Reagent 3) to the section [select the appropriate secondary antibody according to the primary antibody species] and incubate at room temperature for 30-60 minutes.
9. Wash: 3 × 3–5 min with PBST.
10. Chromogenic reaction: Apply DAB substrate (diluted 1:50, e.g., 1 µL DAB concentrate—briefly centrifuge Reagent 5 before use—mixed with 49 µL DAB diluent, Reagent 6) onto the tissue section. Prepare the working solution immediately before use and consume within 30 min. Monitor color development under a light microscope and terminate the reaction promptly once optimal signal intensity is achieved.
11. Wash: 3 × 3–5 min PBST.
12. Blocking: Repeat blocking step with goat serum (Reagent 2) for 30 min at 37 °C to saturate free binding sites.
13. Primary antibody incubation: After determining the optimal dilution, evenly overlay the tissue section with the working antibody solution and incubate overnight at 4 °C or for 1 h at 37 °C.
14. Wash: 3 × 3–5 min PBST.
15. Secondary antibody incubation: Add Polymer-HRP-labeled goat anti-rabbit IgG (Reagent 4) to the section [select the appropriate secondary antibody according to the primary antibody species] and incubate at room temperature for 30-60 minutes.
16. Wash: 3 × 3–5 min PBST.

IHC Double staining kit(Polymer-HRP)(DAB/AEC, For Mouse/Rabbit Primary Antibodies)

Catalog No.: PA10002

-
17. Chromogenic reaction: Apply AEC substrate (prepared fresh at 1:20; e.g., centrifuge Reagents 7 and 8 briefly, then combine 10 μ L AEC chromogen A, 10 μ L AEC chromogen B, and 180 μ L AEC chromogen C) onto the tissue section. Use the working solution within 30 min of preparation. Monitor color development under a light microscope and terminate the reaction promptly once the desired signal intensity is reached.
 18. Counterstain: Stain nuclei with Mayer' s hematoxylin (aqueous formulation; avoid alcohol-containing solutions).
 19. Mounting: Use an aqueous mounting medium (e.g., glycerol gelatin). Do not dehydrate through ethanol or clear in xylene.
 20. Microscopy: Visualize dual-color staining by bright-field microscopy.

Notes

- 1.Store all reagents at 2–8 °C; do not freeze.
- 2.Optimize primary antibody dilutions and incubation times empirically.
- 3.AEC precipitate is alcohol-soluble; avoid alcoholic counterstains, differentiation solutions, ethanol dehydration, and xylene clearing. Use only aqueous mountants.
- 4.Protect AEC substrate and final slides from light to prevent chromogen fading.
- 5.Wear appropriate personal protective equipment (lab coat, gloves, mask).
- 6.For research use only; not intended for diagnostic or therapeutic applications.

This product is for research use only!