

## Mouse CD4<sup>+</sup> T Cell Isolation Kit (negative isolation)

Catalog No.: BRC0001

### Product Description

The Mouse CD4<sup>+</sup> T Cell Isolation Kit is designed to isolate CD4<sup>+</sup> T cells from mouse splenocytes or single-cell suspensions of other tissues via negative selection. The principle involves labeling non-target cells (non- CD4<sup>+</sup> T cells) with biotin-conjugated monoclonal antibodies, followed by depletion using streptavidin-conjugated magnetic beads. This process enriches mouse CD4<sup>+</sup> T cells. A magnetic separator is required for the procedure.

### Kit Components

No.	Reagent Name	Size1 (For 10 <sup>9</sup> cell)	Size2 (For 5×10 <sup>8</sup> cell)
Reagent1	Biotin-Antibody Mix	200µl	100µl
Reagent2	EnkiBeads™ Streptavidin	2 mL	1ml

### Storage Conditions and Shelf Life

**Storage Conditions:** Store at 2–8°C. DO NOT FREEZE.

**Shelf Life:** 2 years

### Applications

This kit is suitable for isolating CD4<sup>+</sup> T cells from mouse spleen and lymph nodes.

### Protocol

#### Example: Isolation of CD4<sup>+</sup> T Cells from Mouse Spleen

1. Prepare single-cell suspension: Mince spleen through a 70 µm cell strainer. Rinse the strainer with ice-cold PBS and collect the cell suspension in a 50 mL centrifuge tube. Centrifuge at 500 g for 5 min.
2. Remove supernatant after centrifugation, add 5 mL erythrocyte lysis buffer (Cat. No. RC0011), and lyse at room temperature for 5 min. Add 20 mL PBS and centrifuge at 500 g for 5 min.

**Note: Erythrocyte lysis conditions may vary depending on the specific lysis buffer used. Trace residual erythrocytes will not affect subsequent isolation or cell purity.**

3. After centrifugation, remove supernatant, resuspend splenocytes in PBS, filter through a 70 µm cell strainer, and count cells. Centrifuge at 500 g for 5 min.

**Note: The cell suspension must be filtered through a cell strainer to remove tissue debris and cell aggregates; otherwise, isolation purity will be compromised.**

4. After centrifugation, remove supernatant and resuspend cells in isolation buffer at a density of 1×10<sup>8</sup> cells/mL.

**Note: Isolation buffer is PBS containing 2 mM EDTA and 2% fetal bovine serum (FBS), or PBS containing 2 mM EDTA and 0.5% BSA. Filter-sterilize through a 0.22 µm membrane before use.**

5. Transfer 100 µL of cell suspension (1×10<sup>7</sup> cells) to the bottom of a sterile flow cytometry tube,

add 2  $\mu\text{L}$  Biotin-Antibody Mix, mix well, and incubate at 4°C for 10 min.

**Note: When adding cell suspension, pipette directly to the bottom of the tube; avoid dispensing along the tube wall. Depending on the magnetic separator used, centrifuge tubes may also be used for cell isolation. For larger cell numbers, increase Biotin-Antibody Mix proportionally.**


6. After incubation, add 20  $\mu\text{L}$  washed EnkiBeads™ Streptavidin to the tube (magnetic beads must be washed before use: vortex to resuspend beads, pipette the required volume into a 1.5 mL centrifuge tube, add 1 mL isolation buffer, centrifuge at 10,000 g for 1 min, and remove supernatant. Repeat wash once with 1 mL isolation buffer, then resuspend beads in the same original volume of isolation buffer. For example, if 20  $\mu\text{L}$  of beads were taken for washing, resuspend in 20  $\mu\text{L}$  isolation buffer after washing), mix well, and incubate at 4°C for 10 min.

**Note: For larger cell numbers, increase EnkiBeads™ Streptavidin proportionally. For example, to isolate  $5 \times 10^7$  cells, add 10  $\mu\text{L}$  Biotin-Antibody Mix and 100  $\mu\text{L}$  EnkiBeads™ Streptavidin to 500  $\mu\text{L}$  cell suspension. For fewer than  $1 \times 10^7$  cells, adjust cell suspension volume to 100  $\mu\text{L}$  and add 2  $\mu\text{L}$  Biotin-Antibody Mix and 20  $\mu\text{L}$  EnkiBeads™ Streptavidin.**

7. After incubation, add 2.5 mL isolation buffer to the tube and mix gently by pipetting up and down 5 times (avoid vigorous shaking or inversion).
8. Place the tube containing cells on the magnetic separator and let stand for 5 min.
9. Gently pour the cell suspension into a sterile centrifuge tube (do not remove the tube from the magnetic separator during pouring). This suspension contains the purified mouse CD4<sup>+</sup> T cells. Centrifuge at 500 g for 5 min, then remove supernatant and collect cells.
10. Wash cells as required for your experiment, then resuspend in appropriate buffer or culture medium for subsequent molecular or cellular biology applications.

## Isolation Performance

CD4<sup>+</sup> T cells were isolated from C57BL/6 mouse splenocytes. Pre- and post-isolation cells were stained with FITC anti-mouse CD4 antibody (clone GK1.5) and analyzed by flow cytometry. CD4<sup>+</sup> T cell purity increased from 20.3% to 98.8%.



## Important Notes

1. Avoid freezing magnetic beads and antibody mix during storage and use.
2. Use low-retention pipette tips and centrifuge tubes to minimize loss of beads and antibodies due to adsorption.
3. This product should be used with a magnetic separator.
4. For research use only.