

Tris-MOPS-SDS

Catalog #: RA10025

Product Overview

This product is white granules with main components of Tris, MOPS, EDTA and SDS. Each sachet of granules can be used to prepare 500mL of electrophoresis running buffer.

Product Information

Product Name	Storage conditions	Specification
Tris-MOPS-SDS	RT	500mL×20
Manual	-	1pcs

Store at 2-8°C ; shelf life is 3 year.

Product Advantages

- **Universal type:** Gel preparation kits and pre-made gels that are compatible with various systems .
- **Highly efficient and time-saving:** 15-35 min rapid electrophoresis.
- **Recyclable:** Can be reused 2-3 times.

Operating procedures

1. Pour 500mL of ddH₂O into a clean beaker;
2. Add 1 sachet of electrophoresis buffer granules to the deionized water and stir until completely dissolved, which is the ready-to-use electrophoresis running buffer .

Precautions

1. This product is for research use only.