

Product Name: cIAP2 (1F12) Mouse Monoclonal Antibody
Catalog #: AMM03485

Summary

Production Name	cIAP2 (1F12) Mouse Monoclonal Antibody
Description	Primary antibody
Host	Mouse
Application	WB
Reactivity	Human,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG1
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	BIRC3
Alternative Names	cIAP2; AIP1; API2; BIRC3; CIAP2; HAIP1; IAP homolog C; RNF49; MIHC; MALT2; Apoptosis inhibitor 2; IAP2
Gene ID	330
SwissProt ID	Q13489

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 68 kDa; Observed MW: 72 kDa

Product Name: cIAP2 (1F12) Mouse Monoclonal Antibody
Catalog #: AMM03485


Background

Apoptotic suppressor. The BIR motifs region interacts with TNF receptor associated factors 1 and 2 (TRAF1 and TRAF2) to form an heteromeric complex, which is then recruited to the tumor necrosis factor receptor 2 (TNFR2).


Research Area

Cell Biology

Image Data



Western blot analysis of cIAP2 in Ramos, COS7, Raji and Daudi lysates using cIAP2 antibody.



Western blot analysis of cIAP2 (1F12) in Ramos, COS7 and Raji lysates using cIAP2 antibody.

Note

For research use only.