

Product Name: Neurofascin Rabbit Monoclonal Antibody
Catalog #: AMRe01464

Summary

Production Name	Neurofascin Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	Nfasc
Alternative Names	Neurofascin
Gene ID	269116.0
SwissProt ID	Q810U3

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 138 kDa; Observed MW: 186,155,75 kDa

Product Name: Neurofascin Rabbit Monoclonal Antibody
Catalog #: AMRe01464


Background

Cell adhesion, ankyrin-binding protein which may be involved in neurite extension, axonal guidance, synaptogenesis, myelination and neuron-glial cell interactions.


Research Area

Neuroscience

Image Data



Western blot analysis of Neurofascin in Human placenta lysates using Neurofascin antibody.



Western blot analysis of Neurofascin in mouse Brain, rat Brain, CSO7, HeLa lysates using Neurofascin antibody.

Note

For research use only.