

Summary

Production Name	ATG3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	ATG3
Alternative Names	APG3; APG3L; PC3-96; APG3-LIKE
Gene ID	64422
SwissProt ID	Q9NT62

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
Molecular Weight	Calculated MW: 36 kDa; Observed MW: 40 kDa

Background

Product Name: ATG3 Rabbit Monoclonal Antibody
Catalog #: AMRe01690




E2-like enzyme involved in autophagy and mitochondrial homeostasis. Catalyzes the conjugation of ATG8-like proteins (GABARAP, GABARAPL1, GABARAPL2 or MAP1LC3A) to phosphatidylethanolamine (PE). Catalyzing the conjugation of ATG12 to itself, ATG12 conjugation to ATG3 playing a role in mitochondrial homeostasis but not in autophagy. ATG7 (E1-like enzyme) facilitates this reaction by forming an E1-E2 complex with ATG3.


Research Area

Cell Biology


Image Data



Western blot analysis of ATG3 in K562, rat Brain, C6, 3T3, HeLa lysates using ATG3 antibody.



Immunocytochemistry analysis of ATG3(green) in HeLa using ATG3 antibody, and DAPI(blue)



Product Name: ATG3 Rabbit Monoclonal Antibody
Catalog #: AMRe01690



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ATG3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.