

Product Name: Bora Rabbit Monoclonal Antibody
Catalog #: AMRe01734



Summary

Production Name	Bora Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	BORA
Alternative Names	C13orf34
Gene ID	79866
SwissProt ID	Q6PGQ7

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 61 kDa; Observed MW: 61 kDa

Background

Product Name: Bora Rabbit Monoclonal Antibody
Catalog #: AMRe01734




Required for the activation of AURKA at the onset of mitosis.


Research Area

Cell Biology

Image Data



Western blot analysis of Bora in mouse liver lysates using Bora antibody.



Western blot analysis of Bora in K562, C6 lysates using Bora antibody.

Note

For research use only.