

Product Name: Calretinin Rabbit Monoclonal Antibody
Catalog #: AMRe01750



Summary

Production Name	Calretinin Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CALB2
Alternative Names	Calretinin; CR; 29 kDa calbindin
Gene ID	794
SwissProt ID	P22676

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200
Molecular Weight	Calculated MW: 32 kDa; Observed MW: 32 kDa

Background

Product Name: Calretinin Rabbit Monoclonal Antibody
Catalog #: AMRe01750




Calretinin is a calcium-binding protein which is abundant in auditory neurons.


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of Calretinin in rat Brain lysates using Calretinin antibody.



Immunocytochemistry analysis of Calretinin(green) in HeLa using Calretinin antibody, and DAPI(blue)

Note

For research use only.