

Product Name: Caspase 7 Rabbit Monoclonal Antibody
Catalog #: AMRe01764



Summary

Production Name	Caspase 7 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF,IP
Reactivity	Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	Casp7
Alternative Names	Mch3; CMH-1; mCASP-7; AI314680; ICE-IAP3; caspase-7
Gene ID	12369.0
SwissProt ID	P97864

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 34 kDa

Background

Product Name: Caspase 7 Rabbit Monoclonal Antibody
Catalog #: AMRe01764




Involved in the activation cascade of caspases responsible for apoptosis execution. Cleaves and activates sterol regulatory element binding proteins (SREBPs). Overexpression promotes programmed cell death .

Research Area

Cell Biology

Image Data



Western blot analysis of Caspase7 in C6, 3T3 lysates using Caspase 7 antibody.

Note

For research use only.