

Product Name: CLPTM1 Rabbit Monoclonal Antibody
Catalog #: AMRe01835



Summary

Production Name	CLPTM1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CLPTM1
Alternative Names	Cleft lip and palate transmembrane protein 1
Gene ID	1209
SwissProt ID	O96005

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 76 kDa; Observed MW: 90 kDa

Background

Product Name: CLPTM1 Rabbit Monoclonal Antibody
Catalog #: AMRe01835




May play a role in T-cell development.

Research Area

Cell Biology

Image Data



Western blot analysis of CLPTM1 in 293T, U87-MG, Jurkat lysates using CLPTM1 antibody.

Note

For research use only.