

Product Name: DEK Rabbit Monoclonal Antibody
Catalog #: AMRe01902



Summary

Production Name	DEK Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	DEK
Alternative Names	Dek; DEK oncogene; Protein DEK
Gene ID	7913
SwissProt ID	P35659

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 50 kDa

Background

Product Name: DEK Rabbit Monoclonal Antibody
Catalog #: AMRe01902




Involved in chromatin organization.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of DEK in K562, MCF-7, HeLa lysates using DEK antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using DEK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.