

Product Name: DHRS1 Rabbit Monoclonal Antibody
Catalog #: AMRe01909



Summary

Production Name	DHRS1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	DHRS1
Alternative Names	dehydrogenase/reductase 1; SDR19C1
Gene ID	115817
SwissProt ID	Q96LJ7

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 34 kDa

Background

Product Name: DHRS1 Rabbit Monoclonal Antibody
Catalog #: AMRe01909




NADPH-dependent oxidoreductase which catalyzes the reduction of steroids (estrone, androstene-3,17-dione and cortisone) as well as prostaglandin E1, isatin and xenobiotics in vitro.

Research Area

Cell Biology

Image Data



Western blot analysis of DHRS1 in C6 lysates using DHRS1 antibody.

Note

For research use only.