

Product Name: DNAJC8 Rabbit Monoclonal Antibody
Catalog #: AMRe01913



Summary

Production Name	DNAJC8 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,IP
Reactivity	Human,Mouse,Rat,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	DNAJC8
Alternative Names	SPF31; HSPC331
Gene ID	22826
SwissProt ID	O75937

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 30 kDa; Observed MW: 30 kDa

Background

Product Name: DNAJC8 Rabbit Monoclonal Antibody
Catalog #: AMRe01913




Suppresses polyglutamine (polyQ) aggregation of ATXN3 in neuronal cells (PubMed:27133716).


Research Area

Signal Transduction

Image Data



Western blot analysis of DjC8 in mouse brain, rat lung lysates using DNAJC8 antibody.



Western blot analysis of DjC8 in Jurkat, HeLa, COS7, Raji lysates using DjC8 antibody.

Note

For research use only.