

Product Name: FAM50A Rabbit Monoclonal Antibody
Catalog #: AMRe01972



Summary

Production Name	FAM50A Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,IP
Reactivity	Human,Mouse,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	FAM50A
Alternative Names	family with sequence similarity 50 member A; 9F; XAP5; HXC26; HXC-26; DXS9928E
Gene ID	9130
SwissProt ID	Q14320

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 40 kDa

Background

Product Name: FAM50A Rabbit Monoclonal Antibody
Catalog #: AMRe01972




May be a DNA-binding protein or transcriptional factor.


Research Area

Cell Biology


Image Data



Western blot analysis of FAM50A/XAP in HeLa, CHO-K1, C6, Jurkat lysates using FAM50A antibody.




Western blot analysis of FAM50A/XAP in Smooth Muscle, HCT116, HeLa lysates using FAM50A/XAP antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using FAM50A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: FAM50A Rabbit Monoclonal Antibody
Catalog #: AMRe01972



Immunocytochemistry analysis of FAM50A (green) in HeLa using FAM50A antibody, and DAPI (blue).

Note

For research use only.