

Product Name: FANCB Rabbit Monoclonal Antibody
Catalog #: AMRe01974



Summary

Production Name	FANCB Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	FANCB
Alternative Names	Fanconi anemia complementation group B; FA2; FAB; FACB; FAAP90; FAAP95
Gene ID	2187
SwissProt ID	Q8NB91

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 98 kDa; Observed MW: 98 kDa

Background

Product Name: FANCB Rabbit Monoclonal Antibody
Catalog #: AMRe01974




DNA repair protein required for FANCD2 ubiquitination.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of FANCB in 3T3, HeLa lysates using FANCB antibody.

Note

For research use only.