

**Product Name: Glucose Transporter GLUT3 Rabbit
Monoclonal Antibody**
Catalog #: AMRe02039



Summary

Production Name	Glucose Transporter GLUT3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	SLC2A3
Alternative Names	SLC2A3; GLUT3; Solute carrier family 2; facilitated glucose transporter member 3; Glucose transporter type 3; brain; GLUT-3
Gene ID	6515
SwissProt ID	P11169

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 54 kDa

**Product Name: Glucose Transporter GLUT3 Rabbit
Monoclonal Antibody**
Catalog #: AMRe02039




Background

Facilitative glucose transporter. Probably a neuronal glucose transporter.


Research Area

Signal Transduction

Image Data



Western blot analysis of Glucose Transporter GLUT3 in mouse brain, mouse heart lysates using Glucose Transporter GLUT3 antibody.



Western blot analysis of Glucose Transporter GLUT3 in Hela, C6, rat Brain lysates using Glucose Transporter GLUT3 antibody.

Note

For research use only.