

**Product Name: Glutathione S Transferase mu 1 Rabbit
Monoclonal Antibody
Catalog #: AMRe02045**



Summary

Production Name	Glutathione S Transferase mu 1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	GSTM1
Alternative Names	MU; H-B; GST1; GTH4; GTM1; MU-1; GSTM1-1; MGC26563; GSTM1a-1a; GSTM1b-1b
Gene ID	2944
SwissProt ID	P09488

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 26 kDa; Observed MW: 26 kDa

**Product Name: Glutathione S Transferase mu 1 Rabbit
Monoclonal Antibody
Catalog #: AMRe02045**




Background

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Research Area

Signal Transduction

Image Data



Western blot analysis of Glutathione S Transferase mu 1 in Jurkat, CHO-K1, HeLa lysates using Glutathione S Transferase mu 1 antibody.

Note

For research use only.