

Product Name: GPRC5B Rabbit Monoclonal Antibody
Catalog #: AMRe02061



Summary

Production Name	GPRC5B Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	GPRC5B
Alternative Names	RAIG2; RAIG-2
Gene ID	51704
SwissProt ID	Q9NZH0

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 45 kDa; Observed MW: 45 kDa

Background

Product Name: GPRC5B Rabbit Monoclonal Antibody
Catalog #: AMRe02061




G-protein coupled receptor family C group 5 member B


Research Area

Neuroscience

Image Data



Western blot analysis of GPRC5B in A549 lysates using GPRC5B antibody.



Western blot analysis of GPRC5B in rat Brain, 3T3 lysates using GPRC5B antibody.

Note

For research use only.