

Product Name: hnRNP G T Rabbit Monoclonal Antibody
Catalog #: AMRe02104



Summary

Production Name	hnRNP G T Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	RBMXL2
Alternative Names	HNRNPGT
Gene ID	27288
SwissProt ID	O75526

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 43 kDa

Background

Product Name: hnRNP G T Rabbit Monoclonal Antibody
Catalog #: AMRe02104




RNA-binding motif protein, X-linked-like-2

Research Area

Tags & Cell Markers

Image Data



Western blot analysis of hnRNP GT in K562, rat Brain, C6, 3T3, HeLa lysates using hnRNP G T antibody.

Note

For research use only.