

Product Name: HPDL Rabbit Monoclonal Antibody
Catalog #: AMRe02116



Summary

Production Name	HPDL Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	HPDL
Alternative Names	4-hydroxyphenylpyruvate dioxygenase like; GLOXD1; 4-HPPD-L
Gene ID	84842
SwissProt ID	Q96IR7

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 39 kDa; Observed MW: 39 kDa

Background

Product Name: HPDL Rabbit Monoclonal Antibody
Catalog #: AMRe02116




May have dioxygenase activity.

Research Area

Signal Transduction

Image Data



Western blot analysis of HPDL in HeLa, Caco2, A549 lysates using HPDL antibody.

Note

For research use only.