

Summary

Production Name	KGF Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	FGF7
Alternative Names	Fgf7; Fibroblast growth factor 7; HBGF7; Heparin binding growth factor 7; Keratinocyte growth factor
Gene ID	2252
SwissProt ID	P21781

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 28 kDa

Product Name: KGF Rabbit Monoclonal Antibody
Catalog #: AMRe02192




Background

Growth factor active on keratinocytes. Possible major paracrine effector of normal epithelial cell proliferation.

Research Area

Cardiovascular

Image Data



Western blot analysis of KGF in rat Brain, C6, 3T3, Hela lysates using KGF antibody.

Note

For research use only.