

Product Name: LGALS3BP Rabbit Monoclonal Antibody
Catalog #: AMRe02213



Summary

Production Name	LGALS3BP Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	LGALS3BP
Alternative Names	90K; M2BP; gp90; CyCAP; BTBD17B; MAC-2-BP; TANGO10B
Gene ID	3959
SwissProt ID	Q08380

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 65 kDa; Observed MW: 100 kDa

Background

Product Name: LGALS3BP Rabbit Monoclonal Antibody
Catalog #: AMRe02213




Promotes integrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells.

Research Area

Cancer

Image Data



Western blot analysis of LGALS3BP in HT-1080, HepG2, A549 lysates using LGALS3BP antibody.

Note

For research use only.