

Product Name: Mesothelin Rabbit Monoclonal Antibody
Catalog #: AMRe02247



Summary

Production Name	Mesothelin Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	Msln
Alternative Names	Pre-pro-megakaryocyte-potentiating factor; MPF
Gene ID	56047.0
SwissProt ID	Q61468

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 70 kDa; Observed MW: 40,70 kDa

Background

Product Name: Mesothelin Rabbit Monoclonal Antibody
Catalog #: AMRe02247




Membrane-anchored forms may play a role in cellular adhesion.

Research Area

Tags & Cell Markers

Image Data



Western blot analysis of Mesothelin in 3T3, mouse lung lysates using Mesothelin antibody.

Note

For research use only.