

Product Name: MKLP1 Rabbit Monoclonal Antibody
Catalog #: AMRe02259



Summary

Production Name	MKLP1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	KIF23
Alternative Names	CHO1; KNSL5; MKLP1; MKLP-1
Gene ID	9493
SwissProt ID	Q02241

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 110 kDa; Observed MW: 110 kDa

Background

Product Name: MKLP1 Rabbit Monoclonal Antibody
Catalog #: AMRe02259




Plus-end-directed motor enzyme that moves antiparallel microtubules in vitro. Localizes to the interzone of mitotic spindles.

Research Area

Cell Biology

Image Data



Western blot analysis of MKLP1 in HeLa lysates using MKLP1 antibody.

Note

For research use only.