

Product Name: MS4A14 Rabbit Monoclonal Antibody
Catalog #: AMRe02275



Summary

Production Name	MS4A14 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	MS4A14
Alternative Names	Testis development protein NYD-SP21
Gene ID	84689
SwissProt ID	Q96JA4

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 77 kDa; Observed MW: 77 kDa

Background

Product Name: MS4A14 Rabbit Monoclonal Antibody
Catalog #: AMRe02275




Membrane-spanning 4-domains subfamily A member 14

Research Area

Signal Transduction

Image Data



Western blot analysis of MS4A14 in Jurkat, rat Brain, C6 lysates using MS4A14 antibody.

Note

For research use only.