

**Product Name: MTH1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe02282**



## Summary

<b>Production Name</b>	MTH1 Rabbit Monoclonal Antibody
<b>Description</b>	Recombinant Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Human

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	NUDT1
<b>Alternative Names</b>	MTH1; NUDT 1
<b>Gene ID</b>	4521
<b>SwissProt ID</b>	P36639

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000
<b>Molecular Weight</b>	Calculated MW: 18 kDa; Observed MW: 18 kDa

## Background

**Product Name: MTH1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe02282**




Antimutagenic. Acts as a sanitizing enzyme for oxidized nucleotide pools, thus suppressing cell dysfunction and death induced by oxidative stress. Hydrolyzes 8-oxo-dGTP, 8-oxo-dATP and 2-OH-dATP, thus preventing misincorporation of oxidized purine nucleoside triphosphates into DNA and subsequently preventing A:T to C:G and G:C to T:A transversions.

## Research Area

Epigenetics and Nuclear Signaling

## Image Data



Western blot analysis of MTH1 in K562 lysates using MTH1 antibody.

## Note

For research use only.