

Product Name: MTX3 Rabbit Monoclonal Antibody
Catalog #: AMRe02287



Summary

Production Name	MTX3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	MTX3
Alternative Names	Metaxin-3
Gene ID	345778
SwissProt ID	Q5HYI7

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 35 kDa

Background

Product Name: MTX3 Rabbit Monoclonal Antibody
Catalog #: AMRe02287




Could function in transport of proteins into the mitochondrion.


Research Area

Signal Transduction

Image Data



Western blot analysis of MTX3 in rat heart lysates using MTX3 antibody.



Western blot analysis of MTX3 in mouse brain lysates using MTX3 antibody.

Note

For research use only.