

Product Name: Mucin 4 Rabbit Monoclonal Antibody
Catalog #: AMRe02289



Summary

Production Name	Mucin 4 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	MUC4
Alternative Names	Mucin-4; Ascites sialoglycoprotein; ASGP; Pancreatic adenocarcinoma mucin; Testis mucin; Tracheobronchial mucin
Gene ID	4585
SwissProt ID	Q99102

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 542 kDa; Observed MW: 120 kDa

Product Name: Mucin 4 Rabbit Monoclonal Antibody
Catalog #: AMRe02289




Background

May play a role in tumor progression. Ability to promote tumor growth may be mainly due to repression of apoptosis as opposed to proliferation.

Research Area

Microbiology

Image Data



Western blot analysis of MUC4 in PC-3 lysates using Mucin 4 antibody.

Note

For research use only.