

Product Name: PCSK9 Rabbit Monoclonal Antibody
Catalog #: AMRe02409



Summary

Production Name	PCSK9 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PCSK9
Alternative Names	HCHOLA3; LDLQC1; NARC1; Neural apoptosis regulated convertase 1; PC9; PCSK9; Proprotein convertase 9
Gene ID	255738
SwissProt ID	Q8NBP7

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 74 kDa; Observed MW: 65,80 kDa

Product Name: PCSK9 Rabbit Monoclonal Antibody
Catalog #: AMRe02409




Background

May be implicated in the differentiation of cortical neurons and may play a role in cholesterol homeostasis.

Research Area

Signal Transduction

Image Data



Western blot analysis of PCSK9 in HeLa lysates using PCSK9 antibody.

Note

For research use only.