

Product Name: RPL10A Rabbit Monoclonal Antibody
Catalog #: AMRe02552



Summary

Production Name	RPL10A Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	RPL10A
Alternative Names	RPL10A; NEDD6; 60S ribosomal protein L10a; CSA-19; Neural precursor cell expressed developmentally down-regulated protein 6; NEDD-6
Gene ID	4736
SwissProt ID	P62906

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 25 kDa; Observed MW: 25 kDa

Product Name: RPL10A Rabbit Monoclonal Antibody
Catalog #: AMRe02552




Background

Component of the large ribosomal subunit.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of RPL10A in K562, C6, 3T3, HeLa lysates using RPL10A antibody.

Note

For research use only.