

Product Name: Striatin 4 Rabbit Monoclonal Antibody
Catalog #: AMRe02645



Summary

Production Name	Striatin 4 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	STRN4
Alternative Names	ZIN; PPP2R6C; zinedin
Gene ID	29888
SwissProt ID	Q9NRL3

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 81 kDa; Observed MW: 81 kDa

Background

Product Name: Striatin 4 Rabbit Monoclonal Antibody
Catalog #: AMRe02645




Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein. Miscellaneous The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).

Research Area

Neuroscience

Image Data



Western blot analysis of Striatin 4 in HeLa lysates using Striatin 4 antibody.

Note

For research use only.