

Product Name: TFG Rabbit Monoclonal Antibody
Catalog #: AMRe02721



Summary

Production Name	TFG Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	TFG
Alternative Names	TF6; HMSNP; SPG57; TRKT3
Gene ID	10342
SwissProt ID	Q92734

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 55 kDa

Background

Product Name: TFG Rabbit Monoclonal Antibody
Catalog #: AMRe02721




Plays a role in the normal dynamic function of the endoplasmic reticulum (ER) and its associated microtubules.


Research Area

Immunology

Image Data



Western blot analysis of TRK fused gene in K562, rat Brain, C6, 3T3, Hela lysates using TFG antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TRK fused gene antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.