

Product Name: UFC1 Rabbit Monoclonal Antibody
Catalog #: AMRe02742



Summary

Production Name	UFC1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	UFC1
Alternative Names	NEDSG; HSPC155
Gene ID	51506
SwissProt ID	Q9Y3C8

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	Calculated MW: 19 kDa; Observed MW: 19 kDa

Background

Product Name: UFC1 Rabbit Monoclonal Antibody
Catalog #: AMRe02742




E2-like enzyme which forms an intermediate with UFM1 via a thioester linkage.


Research Area

Cell Biology

Image Data



Western blot analysis of UFC1 in K562, rat Brain, C6, 3T3 lysates using UFC1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using UFC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.