

Summary

Production Name	USP7 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	USP7 Deubiquitinating enzyme 7; HAUSP; herpes virus associated; Herpes virus associated ubiquitin specific protease; Herpesvirus-associated ubiquitin-specific protease; TEF 1; tef-1; TEF1; Ubiquitin carboxyl terminal hydrolase 7; Ubiquitin carboxyl-terminal hydrolase 7; Ubiquitin specific peptidase 7 (herpes virus associated); Ubiquitin specific peptidase 7; Ubiquitin specific peptidase 7 herpes virus associated; Ubiquitin specific processing protease 7; Ubiquitin specific protease 7 (herpes virus associated); Ubiquitin Specific Protease 7; Ubiquitin specific protease 7 herpes virus associated; Ubiquitin thioesterase 7; ubiquitin thiolesterase 7; Ubiquitin-specific-processing protease 7; UBP 7; UBP-7; UBP7; UBP7_HUMAN; USP 7; usp-7; USP7; VMW110-ASSOCIATED PROTEIN; 135-KD.
Alternative Names	
Gene ID	7874

Product Name: USP7 Rabbit Monoclonal Antibody
Catalog #: AMRe02750



SwissProt ID

Q93009

Application

Dilution Ratio

WB: 1/500-1/1000 IHC: 1/50-1/100

Molecular Weight

Calculated MW: 128 kDa; Observed MW: 128 kDa


Background

Hydrolase that deubiquitinates target proteins such as FOXO4, p53/TP53, MDM2, ERCC6, DNMT1, UHRF1, PTEN and DAXX (PubMed:11923872, PubMed:15053880, PubMed:16964248, PubMed:18716620, PubMed:25283148).


Research Area

Cell Biology

Image Data



Western blot analysis of USP7 in C6 lysates using USP7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HAUSP/USP7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: USP7 Rabbit Monoclonal Antibody
Catalog #: AMRe02750



Note

For research use only.