

## Summary

<b>Production Name</b>	VAMP2 Rabbit Monoclonal Antibody
<b>Description</b>	Recombinant Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ICC/IF,IP
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	VAMP2
<b>Alternative Names</b>	vesicle-associated membrane protein 2 (synaptobrevin 2); SYB2; VAMP-2
<b>Gene ID</b>	6844
<b>SwissProt ID</b>	P63027

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
<b>Molecular Weight</b>	Calculated MW: 13 kDa; Observed MW: 13 kDa

## Background

**Product Name: VAMP2 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe02751**




The SNARE complex is involved in vesicular transport and membrane fusion, a process regulated by calcium. In neurons, VAMP2 is predominantly inserted in presynaptic vesicle membranes. Assembly of VAMP2 with the plasma membrane SNAREs syntaxin 1 and SNAP25 is a key event necessary for membrane fusion and neurotransmitter release.


## Research Area

Signal Transduction


## Image Data



Immunocytochemistry analysis of VAMP2 (green) in Jurkat using VAMP2 antibody, and DAPI (blue).



Western blot analysis of VAMP2 in mouse brain lysates using VAMP2 antibody.



Western blot analysis of VAMP2 in rat Brain lysates using VAMP2 antibody

**Product Name: VAMP2 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe02751**

---



**Note**

For research use only.