

Summary

Production Name	CAP2 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CAP2
Alternative Names	CAP 2; Adenylyl cyclase-associated protein 2
Gene ID	10486
SwissProt ID	P40123

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Background

Product Name: CAP2 Rabbit Monoclonal Antibody
Catalog #: AMRe02821




May have a regulatory bifunctional role.

Research Area

Signal Transduction

Image Data



Western blot analysis of CAP2 in HeLa, CHO-K1 lysates using CAP2 antibody.

Note

For research use only.