
Summary

Production Name	RASGRP1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	RASGRP1
Alternative Names	V; RASGRP; hRasGRP1; CALDAG-GEFI; CALDAG-GEFII
Gene ID	10125
SwissProt ID	O95267

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 90 kDa; Observed MW: 90 kDa

Background

Product Name: RASGRP1 Rabbit Monoclonal Antibody
Catalog #: AMRe02827




Functions as a calcium- and diacylglycerol (DAG)-regulated nucleotide exchange factor specifically activating Ras through the exchange of bound GDP for GTP.

Research Area

Neuroscience

Image Data



Western blot analysis of Rasgrp1 in A549, HL-60, U251 lysates using RASGRP1 antibody.

Note

For research use only.