

Product Name: Bcl2 L14 Rabbit Monoclonal Antibody
Catalog #: AMRe02911



Summary

Production Name	Bcl2 L14 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	BCL2L14
Alternative Names	BCLG
Gene ID	79370
SwissProt ID	Q9BZR8

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 37 kDa; Observed MW: 43 kDa

Background

Product Name: Bcl2 L14 Rabbit Monoclonal Antibody
Catalog #: AMRe02911




Plays a role in apoptosis.


Research Area

Cell Biology


Image Data



Western blot analysis of Bcl G/BCL2L14 in C6, 3T3, PC-12, Raw264.7 lysates using Bcl2 L14 antibody.



Western blot analysis of Bcl G/BCL2L14 in HeLa, A549, HL-60, U251, U87-MG lysates using Bcl G/BCL2L14 antibody.



Immunocytochemistry analysis of Bcl G/BCL2L14 (green) in HeLa using Bcl G/BCL2L14 antibody, and DAPI (blue)

Product Name: Bcl2 L14 Rabbit Monoclonal Antibody
Catalog #: AMRe02911



Note

For research use only.