

Product Name: Integrin alpha L Rabbit Monoclonal Antibody
Catalog #: AMRe02932

Summary

Production Name	Integrin alpha L Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	Itgal
Alternative Names	Cd11a; LFA-1; Ly-15; Ly-21; (p180); LFA-1A
Gene ID	16408.0
SwissProt ID	P24063

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	Calculated MW: 128 kDa; Observed MW: 160 kDa

Product Name: Integrin alpha L Rabbit Monoclonal Antibody
Catalog #: AMRe02932




Background

Integrin ITGAL/ITGB2 is a receptor for ICAM1, ICAM2, ICAM3 and ICAM4.


Research Area

Tags & Cell Markers

Image Data



Immunohistochemistry analysis of paraffin-embedded mouse lung using Integrin alpha L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD11a in A549, HL-60 lysates using CD11a antibody.

Note

For research use only.