

Product Name: Cellubrevin Rabbit Monoclonal Antibody
Catalog #: AMRe02941



Summary

Production Name	Cellubrevin Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	VAMP3
Alternative Names	CEB; Cellubrevin; Synaptobrevin 3; VAMP3
Gene ID	9341
SwissProt ID	Q15836

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 11 kDa; Observed MW: 11 kDa

Background

Product Name: Cellubrevin Rabbit Monoclonal Antibody
Catalog #: AMRe02941




SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.


Research Area

Signal Transduction

Image Data



Western blot analysis of Cellubrevin in HeLa, A549, MCF-7, 293T lysates using Cellubrevin antibody.



Immunocytochemistry analysis of Cellubrevin (green) in HeLa using Cellubrevin antibody, and DAPI (blue)

Note

For research use only.