

Product Name: CrkL Rabbit Monoclonal Antibody
Catalog #: AMRe02952



Summary

Production Name	CrkL Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CRKL
Alternative Names	Crk L; Crk like protein; Crkl; Crkol
Gene ID	1399
SwissProt ID	P46109

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 39 kDa

Background

Product Name: CrkL Rabbit Monoclonal Antibody
Catalog #: AMRe02952




May mediate the transduction of intracellular signals.


Research Area

Neuroscience

Image Data



Immunocytochemistry analysis of CrkL (green) in K562 using CrkL antibody, and DAPI (blue).



Western blot analysis of CrkL in Jurkat lysates using CrkL antibody.

Note

For research use only.