

Product Name: DENN Rabbit Monoclonal Antibody
Catalog #: AMRe02968



Summary

Production Name	DENN Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	MADD
Alternative Names	DENN; IG20; RAB3GEP
Gene ID	8567
SwissProt ID	Q8WYG6

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
Molecular Weight	Calculated MW: 183 kDa; Observed MW: 220 kDa

Background

Product Name: DENN Rabbit Monoclonal Antibody
Catalog #: AMRe02968




Guanyl-nucleotide exchange factor that regulates small GTPases of the Rab family.


Research Area

-


Image Data



Western blot analysis of DENN in HeLa, A549, HL-60 lysates using DENN antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using DENN antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of DENN (green) in HeLa using DENN antibody, and DAPI (blue)

Product Name: DENN Rabbit Monoclonal Antibody
Catalog #: AMRe02968



Note

For research use only.