

Product Name: EYA2 Rabbit Monoclonal Antibody
Catalog #: AMRe02991



Summary

Production Name	EYA2 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	EYA2
Alternative Names	EAB1
Gene ID	2139
SwissProt ID	O00167

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 59 kDa; Observed MW: 59 kDa

Background

Product Name: EYA2 Rabbit Monoclonal Antibody
Catalog #: AMRe02991




Functions both as protein phosphatase and as transcriptional coactivator for SIX1, and probably also for SIX2, SIX4 and SIX5.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of EYA2 in HeLa, A549, U2OS lysates using EYA2 antibody.

Note

For research use only.