

Summary

Production Name	FGFR3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	FGFR3
Alternative Names	FGFR3; JTK4; Fibroblast growth factor receptor 3; FGFR-3; CD antigen CD333
Gene ID	2261
SwissProt ID	P22607

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 125 kDa

Background

Product Name: FGFR3 Rabbit Monoclonal Antibody
Catalog #: AMRe02998




Receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.


Research Area

Cardiovascular

Image Data



Western blot analysis of FGFR3 in mouse kidney lysates using FGFR3 antibody.



Western blot analysis of FGFR3 in HeLa, A549, HL-60, U2OS, C6 lysates using FGFR3 antibody.

Note

For research use only.