

Summary

Production Name	LTK Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	LTK
Alternative Names	TYK1
Gene ID	4058
SwissProt ID	P29376

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100
Molecular Weight	Calculated MW: 92 kDa; Observed MW: 92 kDa

Background

Product Name: LTK Rabbit Monoclonal Antibody
Catalog #: AMRe03071




Receptor with a tyrosine-protein kinase activity.


Research Area

-

Image Data



Western blot analysis of LTK in A549, HL-60, U2OS, C6 lysates using LTK antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using LTK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.