

**Product Name: MAPKAP Kinase 3 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe03088**

---



## Summary

<b>Production Name</b>	MAPKAP Kinase 3 Rabbit Monoclonal Antibody
<b>Description</b>	Recombinant Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IP
<b>Reactivity</b>	Human

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	MAPKAPK3
<b>Alternative Names</b>	MAPKAPK3; MAP kinase-activated protein kinase 3; MAPK-activated protein kinase 3; MAPKAP kinase 3; MAPKAP-K3; MAPKAPK-3; MK-3; Chromosome 3p kinase; 3pK
<b>Gene ID</b>	7867
<b>SwissProt ID</b>	Q16644

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IP: 1/20
<b>Molecular Weight</b>	Calculated MW: 43 kDa; Observed MW: 43 kDa

**Product Name: MAPKAP Kinase 3 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe03088**

---




## Background

Stress-activated serine/threonine-protein kinase involved in cytokines production, endocytosis, cell migration, chromatin remodeling and transcriptional regulation.

## Research Area

Signal Transduction

## Image Data



Western blot analysis of MK3 in A549, HL-60 lysates using MAPKAP Kinase 3 antibody.

## Note

For research use only.