

Summary

Production Name	PPCEL Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PREPL
Alternative Names	CMS22; Prolylendopeptidase-like
Gene ID	9581
SwissProt ID	Q4J6C6

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 84 kDa; Observed MW: 72 kDa

Background

Product Name: PPCEL Rabbit Monoclonal Antibody
Catalog #: AMRe03159




Serine peptidase whose precise substrate specificity remains unclear.


Research Area

-


Image Data



Immunocytochemistry analysis of PPCEL (green) in U87-MG using PPCEL antibody, and DAPI (blue).



Western blot analysis of PPCEL/PREPL in mouse brain lysates using PPCEL/PREPL antibody.



Western blot analysis of PPCEL/PREPL in HeLa, C6, Ramos lysates using PPCEL/PREPL antibody

Product Name: PPCEL Rabbit Monoclonal Antibody
Catalog #: AMRe03159



Note

For research use only.