

Product Name: Protein CASP Rabbit Monoclonal Antibody
Catalog #: AMRe03165

Summary

Production Name	Protein CASP Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CUX1
Alternative Names	CDP; CUX; p75; CASP; CDP1; COY1; Clox; GDDI; p100; p110; p200; CUTL1; GOLIM6; CDP/Cut; Cux/CDP; Nbla10317
Gene ID	1523
SwissProt ID	Q13948

Application

Dilution Ratio	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Molecular Weight	Calculated MW: 77 kDa; Observed MW: 77 kDa

Product Name: Protein CASP Rabbit Monoclonal Antibody
Catalog #: AMRe03165


Background

May be involved in intra-Golgi retrograde transport.


Research Area

-


Image Data



Immunocytochemistry analysis of Protein CASP (green) in A549 using Protein CASP antibody, and DAPI (blue).



Western blot analysis of Protein CASP in 3T3, Raw264.7 lysates using Protein CASP antibody.





Product Name: Protein CASP Rabbit Monoclonal Antibody
Catalog #: AMRe03165

Western blot analysis of Protein CASP in Hela, A549, HL-60, U2OS, C6 lysates using Protein CASP antibody

Note

For research use only.