

Product Name: SEC24D Rabbit Monoclonal Antibody
Catalog #: AMRe03197



Summary

Production Name	SEC24D Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	SEC24D
Alternative Names	CLCRP2
Gene ID	9871
SwissProt ID	O94855

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 113 kDa; Observed MW: 113 kDa

Background

Product Name: SEC24D Rabbit Monoclonal Antibody
Catalog #: AMRe03197




Component of the coat protein complex II (COPII) which promotes the formation of transport vesicles from the endoplasmic reticulum (ER).


Research Area

-

Image Data



Western blot analysis of SEC24D in Raw264.7, C6 lysates using SEC24D antibody.



Western blot analysis of SEC24D in HeLa, A549, HL-60 lysates using SEC24D antibody.

Note

For research use only.