

Product Name: Serine Racemase Rabbit Monoclonal Antibody
Catalog #: AMRe03200



Summary

Production Name	Serine Racemase Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	SRR
Alternative Names	ILV1; ISO1
Gene ID	63826
SwissProt ID	Q9GZT4

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 37 kDa; Observed MW: 37 kDa

Product Name: Serine Racemase Rabbit Monoclonal Antibody
Catalog #: AMRe03200




Background

Catalyzes the synthesis of D-serine from L-serine.

Research Area

-

Image Data



Western blot analysis of Serine racemase in HeLa, A549, HL-60, U251 lysates using Serine Racemase antibody.

Note

For research use only.